



# Infinity

## All-Purpose Resin-Ionomer Cement

Infinity is a fluoride-releasing, resin-ionomer, biocompatible dual-cure cement. Its high bonding strength, tough fracture resistance, low film thickness and radiopacity make Infinity the ideal material for cementing ALL restorations. Infinity Syringeable is proven to successfully bond to all surfaces—even zirconia.

Infinity is packaged in a dual-barrel syringe with auto-mixing tips. Easy and simple to dispense, Infinity Syringeable requires no hand-mixing or extra equipment, enabling you to produce high quality, dual-cure cementations every time, on every surface.



- Cleans up in seconds
- Self-adhesive
- Self-etch
- Dual-cure formula
- Low film thickness
- High bond strengths to all surfaces

### Applications

- All-porcelain crown and bridge cementation
- PFM and monolithic zirconia restoration cementation
- Amalgam cementation
- Fluoride-releasing base/liner
- Maryland bridge cementation
- Endodontic post-and-core cementations
- Resin or ceramic inlays and onlays

### Infinity Cementation Technique



- 1**  
Clean tooth surfaces  
Prepare Infinity cement for use and dispense into the restoration



- 2**  
Seat the restoration under occlusal forces



- 3**  
Briefly light-cure excess material at the margins for 1-3 seconds until it reaches a gel consistency



- 4**  
Remove excess gelled material with an instrument



- 5**  
Use dental floss to remove excess material interproximally



- 6**  
Continue light-curing process or allow the restoration to set for 2 to 3 minutes

Available through  
 HENRY SCHEIN®  
DENTAL

To Order: 1-800-372-4346 8:00am-8:30pm (et)  
[www.henryscheindental.com](http://www.henryscheindental.com)

Features	Benefits
Auto-mixing syringe	Simple, easy to use, minimal waste, time saving
Insoluble	Prevents marginal wash-out and the need for recementation of crowns
Self-adhesive, self-etch cement	One product for adhesive procedures Reduces inventory, saves money
High bond strengths to all surfaces	Assures cementation success Eliminates need for multiple cement products
Low film thickness	Allows for precise fitting of crowns and inlays/onlays
Contains fluoride	Continuous slow release of fluoride
Dual-cure	Offers complete curing with or without a light
Radiopaque	Highly distinguishable from natural tooth structure in x-rays

### Infinity is perfect for cementing:



PFM



Zirconia



All Ceramic



Inlays & Onlays



Posts



Product	Item #
<b>Infinity Intro Kit</b>	
Infinity Intro Kit	039690110
<i>Contents: 1 – 5 g syringe, auto-mix dispensing tips (3 pc.), instructions</i>	
<b>Infinity Value Kits</b>	
Infinity Syringeable Value Kit	039680100
Infinity SE Syringeable Value Kit	039681200
<i>Contents: 3 – 10 g syringes, auto-mix dispensing tips (60 pc.), instructions</i>	
<b>Infinity Syringeable Kits</b>	
Infinity Syringeable Kit	039670100
Infinity SE Syringeable Kit	039680200
<i>Contents: 1 – 10 g syringe, auto-mix dispensing tips (20 pc.), instructions</i>	
<b>Infinity Jars</b>	
Infinity Original A & B Jars	039620150
<i>Contents: 1 – 8 g A jar, 1 – 8 g B jar, mixing pad, mixing sticks, instructions</i>	
<b>Product</b>	
<b>Ancillary Cement Accessories</b>	
Syringe Mixing Tips 8-Element (10 pc.)	030795525

Available through



To Order: 1-800-372-4346 8:00am–8:30pm (et)  
[www.henryscheindental.com](http://www.henryscheindental.com)