


Get hooked!



Hartzell Instruments make Gutta-Percha removal easy!

Meticulously crafted by Hartzell Instruments, the specialized Gutta-Percha Removal Instruments have a microscopic (.06mm!) retrieval hook or diamond-coated tip to help simplify your endodontic retreatments. Just insert the instrument into the canal, “hook” the Gutta-Percha and remove!

- ✓ Available in Left/Right and Up/Down directional hooks to retrieve Gutta-Percha anywhere in the canal
- ✓ Flexible tips to easily conform to the shape of the canal
- ✓ Available in 18mm and 30mm lengths to accommodate most depth needs
- ✓ Fabricated using highly polished stainless steel for easy sterilization

Hartzell
INSTRUMENTS

Available through

 **HENRY SCHEIN®**
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)
www.henryscheindental.com


DenMat®