



# 50% of patients have gingivitis.

How do you feel about that?

## Feel confident with Hartzell's Gingivitis Therapy Kit

Featuring cone socket instruments!

### Detect. Access. Remove.

Reverse and remedy your patients' gingivitis, while reducing costs and waste of instrument replacement with cone socket tips!



### Kit Includes:



#### SAMJ 1/2 Sickle Scaler

The workhorse of scalers. Accesses and removes all calculus from **supragingival** areas of the teeth.

PN: HZMJ1S & HZMJ2S



#### 1/2 Gracey Anterior Curette

Effectively accesses and removes **subgingival** calculus in all anterior pockets – canine to canine.

PN: HZG1S & HZG2S



#### 13/14 Gracey Posterior Curette

Effectively accesses and removes **subgingival** calculus in all distal posterior pockets – canine to canine.

PN: HZG13S & HZG14S.



#### 15/16 Gracey Posterior Curette

Contemporary design accesses and removes **subgingival** calculus in all mesial, buccal and lingual posterior pockets.

PN: HZG15S & HZG16S.



#### 11/12 Double-Ended Explorer\*

Detects subgingival calculus in all areas of inflammation.

PN: HZXA11/2

\*Double-ended explorer is not available in cone socket tips.



#### Paired with our lightweight, ergonomic "A" handles

Prevents hand, wrist, and elbow fatigue without hindering tactile transfer.

PN: HZCSHA



#### Cone Socket Wrench

Secure each tip into handle with wrench for optimal safety and efficiency.

PN: HZCSWS

Available through

 **HENRY SCHEIN®**  
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)  
[www.henryscheindental.com](http://www.henryscheindental.com)

  
DenMat®