

## Core Paste XP

### Core Build-Up Material

Core Paste XP has the same strength and excellent handling characteristics of original Core Paste, with over 40 years of clinical success. The stackable viscosity allows you to build up without a matrix or core form, and the dual-cure formula gives you the confidence of a self-cure with the option of light-curing to speed the procedure.

Core Paste XP is a creamy, radiopaque paste developed specifically for high strength reliable core build-ups. When there is insufficient tooth structure available for a fixed prosthesis, CorePaste XP provides the added strength and stability needed for long-term clinical success.

- Dual-cure formula
- Stackable viscosity and excellent handling characteristics
- Cuts like dentin
- Radiopaque

### Applications

- Core build-ups
- Integral post and core procedures
- Cementation of posts
- Esthetic posterior restorations



Before



After

Core Paste's stackable viscosity allows build-up without a matrix or core form!

Available through

 **HENRY SCHEIN®**  
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)  
[www.henryscheindental.com](http://www.henryscheindental.com)

Features	Benefits
Excellent handling characteristics	Stackable viscosity, syringes easily and cuts like dentin
Enamel and white shades	Enamel shade enhances the esthetics of translucent restorations White shade is easily distinguished from tooth structure
Fluoride release	Helps to prevent secondary decay
Auto-mixing syringe delivery	Eliminates hand mixing; less waste than cartridge/gun systems
Radiopaque	Easily distinguishable in radiographs

“I’ve been using Core Paste in my practice for many years. The stackable viscosity provides excellent handling properties. Core Paste is easy to use and consistently results in strong and effective restorations.”

— Ian Shuman, DDS, Pasadena, MD

### Core Post Ideal for Core Build-Ups in Single-Tooth Restorations



#### Core Post Value Kit

#033637000

50 posts

Assorted sizes of carbon fiber and glass fiber posts

Product	Item #
<b>Core Paste XP</b>	
Core Paste XP Enamel with Fluoride Intro Kit*	030657110
Core Paste XP Dual-Cure Enamel with Fluoride*	030652110
Core Paste XP Dual-Cure Enamel	030653110
Core Paste XP Dual-Cure White	030654110
Core Paste XP Dual-Cure White with Fluoride*	030655110
<small>Contents: 4 – 10 g syringes (Intro Kit contains 2 syringes only), auto-mixing tips (40 pc.), clear intraoral tips (40 pc.), intraoral tips 18 G (6 pc.), Tenure A &amp; B bonding agent, Etch'N'Seal</small>	

Product	Item #
<b>Core Paste Syringeable</b>	
Core Paste White with Fluoride Self-Cure*	030625100
Core Paste White Self-Cure	030635100
Core Paste White with Fluoride Self-Cure & Gun*	030625500
Core Paste White Self-Cure & Gun	030635500
Core Paste Enamel Dual-Cure	030645100
Core Paste Enamel Dual-Cure & Gun	030645500
<small>Contents: 1 – 50 g cartridge, auto-mixing tips (20 pc.), intraoral tips (20 pc.), mixing pad, mixing tips, instructions</small>	

Product	Item #
<b>Core Paste Jars</b>	
Core Paste White Q/C with Fluoride Kit* – 50 g	030612525
Core Paste White with Fluoride Kit* – 50g	030622100
Core Paste Enamel Kit – 50 g	030616000
Core Paste Enamel with Fluoride Kit* – 50 g	030616025
Core Paste Enamel Q/C Kit – 50 g	030616325
Core Paste Blue S/C Kit – 50 g	030617325
Core Paste White Kit – 50 g	6100
<small>Contents: 2 – 25 g jars Core Paste (A &amp; B), centrix tip/plunger orange/green (25 pc.), mixing pad, mixing sticks, All-Surface Bonding Chart, instructions</small>	

Available through

 **HENRY SCHEIN®**  
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)  
**www.henryscheindental.com**