



Perfectemp II

Temporary Crown & Bridge Material

Perfectemp II is the tried-and-true provisional with a 10-year satisfaction history. This temporary crown and bridge material has the strength to resist breakage and wearing, combined with the shade options and polish-ability to match the esthetics of a natural-looking tooth.

- 6 min set time; 2.5 min removal time
- Easy trimming and polishing for a natural luster
- 5 Vita shades for a perfect match
- Both cartridge and syringe options
- Affordable pricing compared to leading competitors



Features	Benefits
Only 0.45% shrinkage	Perfect-fitting temporaries
Intraoral set time 2:30 and total set time 6:00	Saves time and money
5 shades: A1, A2, A3.5, B1 and B2	Versatility to nearly match every patient
Natural fluorescence	Creates a natural tooth-like appearance
Minimal heat generation	Maintains pulpal integrity
Automix/autocure	Easy to use, no mess
Multiple uses	Single or multiple crowns, bridges, inlays and onlays

Physical Properties	
Compressive strength (MPa)	≥150
Transverse strength (MPa)	≥60
Diametral tensile strength (MPa)	31
Barcol hardness	≥20

Product	Item #
Perfectemp II Syringe – 5 ml	
A1	CR1039
A2	CR1043
A3.5	CR1042
B1	CR1041
BL-2	CR1040
Contents: 2 – 5 ml syringes, mixing tips, instructions	

Product	Item #
Perfectemp II Cartridge – 50 ml	
A1	SPD1300
A2	SPD1301
A3.5	SPD1302
B1	SPD1303
BL-2	SPD1320
Contents: 1 – 50 ml cartridge, mixing tips, instructions	

©2018 Den-Mat Holdings, LLC. All rights reserved. 821266500HS R:0 05/18SN

Available through

To Order: 1-800-372-4346 8:00am–8:30pm (et)
www.henryscheindental.com