

Bite

Bite Registration Material

DenMat offers a complete line of VPS bite registration materials that ensure precision for your lab and meet all of your procedural needs. Easy to extrude and trim, these materials are used for closed, open and double bite registrations. No matter your preference or technique, DenMat has a bite registration material sure to satisfy every clinician.

With Bite you get extreme accuracy and stability, fast set times and easy trimming.



Bite is specifically engineered to consistently deliver dimensionally stable, accurate and rigid bite registrations

- Mousse-like consistency — no resistance to closure, ensuring proper centric occlusion
- Rigid final set — reduces distortion
- Fast set times
- Easy trimming



Available through

 **HENRY SCHEIN[®]**
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)
www.henryscheindental.com



	Vanilla	Chocolate	Clear	MegaBite
Working Time	20 seconds	45 seconds	30 seconds	20 seconds
Total Set Time	55 seconds	1:30 minutes	2:30 minutes	55 seconds
Hardness	shore A 90	shore A 90	shore A 63	shore D > 40
Color	white	brown	clear	teal
Scent	vanilla	chocolate	mint	mint



Vanilla Bite™

- Mousse-like consistency
- Rigid
- Fast set time (0:55)
- Ideal for 1-2 unit cases
- Extremely accurate and dimensionally stable



Chocolate Bite™

- Slower set time (1:30)
- Ideal for taking full-arch quadrant registrations
- Allows for longer working times



Clear Bite™

- Transparent
- Flexible but resilient
- Ideal for temporary matrix (2:30 set time)
- Enables light passage for curing
- Great visibility of voids or bubbles



MegaBite®

- Very rigid VPS bite registration material
- Fast set time (0:55)
- Exceptional accuracy
- Ideal for centric occlusion and centric relation registrations
- Minimized occlusal adjustments

Product	Item #
Bite Registration	
Vanilla Bite	SPD1105
Chocolate Bite	SPD1186
Clear Bite	SPD1104
MegaBite	SPD1370
<i>Contents: 2 – 48 ml cartridges, tips, instructions</i>	

Available through

 **HENRY SCHEIN®**
DENTAL

To Order: 1-800-372-4346 8:00am–8:30pm (et)
www.henryscheindental.com