



SplashMax VPS Impression Material

SplashMax is a VPS impression material that offers instant hydrophilicity, because dead-on accuracy is needed when the material first comes in contact with dentition. In independent testing, SplashMax demonstrated the lowest contact angle almost immediately. It's the ultimate in accuracy, reliability and value!

- Extreme wettability
- Maximum tear strength and elasticity
- Excellent dimensional accuracy
- Four high-contrast viscosities
- Wild berry scent




Available through

 **HENRY SCHEIN®**
DENTAL

To Order: 1-800-372-4346 8:00am-8:30pm (et)
www.henryscheindental.com



Viscosity	Set Time	Color
Extra Lite Body	2:15 & 4:30	Orange
Lite Body	2:15 & 4:30	Green
Medium Body	2:15 & 4:30	Blue
Heavy Body	2:15 & 4:30	Purple



SplashMax features the lowest contact angle within seconds!*



*Data on file.

Product	Item #
SplashMax 2-Pack Cartridges Regular Set (4:30)	
Extra Lite Body	SPD1613
Lite Body	SPD1612
Medium Body	SPD1611
Heavy Body	SPD1610
<i>Contents: 2 – 50 ml cartridges, tips, instructions</i>	
SplashMax 2-Pack Cartridges Half-Time Set (2:15)	
Extra Lite Body	SPD1617
Lite Body	SPD1616
Medium Body	SPD1615
Heavy Body	SPD1614
<i>Contents: 2 – 50 ml cartridges, tips, instructions</i>	

Product	Item #
SplashMax 8-Pack Cartridges Regular Set (4:30)	
Extra Lite Body	SPD1606
Lite Body	SPD1608
<i>Contents: 8 – 50 ml cartridges, instructions</i>	
SplashMax 8-Pack Cartridges Half-Time Set (2:15)	
Extra Lite Body	SPD1607
Lite Body	SPD1609
<i>Contents: 8 – 50 ml cartridges, instructions</i>	

Product	Item #
SplashMax 20-Pack Cartridges Regular Set (4:30)	
Medium Body	SPD1604
Heavy Body	SPD1602
<i>Contents: 20 – 50 ml cartridges, instructions</i>	
SplashMax 20-Pack Cartridges Half-Time Set (2:15)	
Medium Body	SPD1605
Heavy Body	SPD1603
<i>Contents: 20 – 50 ml cartridges, instructions</i>	
SplashMax 2-Pack Jumbo Cartridges Regular Set (4:30)	
Medium Body	SPD1601
Heavy Body	SPD1600
<i>Contents: 2 – 360 ml jumbo cartridges, mixing tips, rotating rings, instructions</i>	

“SplashMax offers hydrophilic and flowable material properties that allow me to capture detailed, bubble-free impressions for all of my restorative cases. The bright colors provide impressions that are easy to read and patients prefer the taste to other materials we have used in the past. SplashMax is the go-to impression material for my practice!”

– Timothy F. Kosinski, DDS, MAGD



Available through



To Order: 1-800-372-4346 8:00am–8:30pm (et)
www.henryscheindental.com